



# Quality Assurance Panel

Jonathan Mosley<sup>1</sup>, Annie Evans<sup>2</sup> and Susan Sumner<sup>3</sup>

<sup>1</sup>National Exposure Research Laboratory, U.S. EPA, Athens, GA, USA

<sup>2</sup>Metabolon, Research Triangle Park, NC, USA

<sup>3</sup>University of North Carolina at Chapel Hill, Chapel Hill, NC, USA

# QA vs QC: American Society for Quality Definitions

- **Quality assurance** consists of that “part of *quality management* focused on providing confidence that *quality requirements* will be fulfilled.” The confidence provided by quality assurance is twofold—internally to management and externally to customers, government agencies, regulators, certifiers, and third parties.
- **Quality control** is that “part of *quality management* focused on fulfilling *quality requirements*.”

<http://asq.org/learn-about-quality/quality-assurance-quality-control/overview/overview.html>

# What is QA??


---

<u>Parameters</u>	<u>Quality Assurance</u>	<u>Quality Control</u>
<b>Definition</b>	QA is a set of activities for ensuring quality in the <b>processes</b> by which products are developed.	QC is a set of activities for ensuring quality in products. The activities focus on identifying defects in the actual products produced.
<b>Focus on</b>	QA aims to prevent defects with a focus on the <b>process</b> used to make the product. It is a <u>proactive</u> quality process.	QC aims to identify (and correct) defects in the finished product. Quality control, therefore, is a reactive process.
<b>Goal</b>	The goal of QA is to improve development and test <b>processes</b> so that defects do not arise when the product is being developed.	The goal of QC is to identify defects after a product is developed and before it's released.

# QA Definition Proposed for Workshop

*QA is the set of procedures that are performed in advance of analysis of samples that are used to improve data quality (e.g. education, analyst competency, method development and documentation)*

# Summary of QA/QC procedures



See A. Evans' **Poster:** Session 1 (Monday 5:15pm - 6:45pm ) and Session 2 (Tuesday 5:15pm - 6:45pm)  
**Poster Group:** Sample Prep and Quality **Poster Number:** 333

# General Ideas

---

- Keep It Simple Silly – if you have nothing in place, start simply
- Part of a quality system is to be able to trace what you are doing... if the system is too complicated or too cumbersome people will not follow it and traceability goes to zero.
- Establish clear criteria
- Check lists are actually pretty useful. They are used in flight cockpits!
- Any person in a group can operate as a QA person at any time. You don't need a dedicated QA person (though this is also an option). "Sure I can review your calculations".

# Pre-Analysis Ideas


- Training Instrument Operators
  - Weighing tests
  - Pipet tests
  - Best practices for notebooks
- Instrument Maintenance and Troubleshooting
  - Tuning and calibration
  - Column Pressure Checks
  - Acceptance criteria for benchmarking
- Checklist – what to do if acceptance criteria are not met
- Standard Operating Procedures (SOP)
  - Maintained in a repository (i.e. needs version control)
  - Cover every step of the analytical process
  - Also for post-analytical processes (e.g. data pre-processing and calculations)

# “Document, document, document ...”

- Notebooks, Log Books (e.g. instrument, sample storage), Temperature Monitoring Sheets, etc.
- List of SOPs (or Method files) that each individual has signed as “read and understood”
- List of SOPs that the trainer has signed off on for each individual.
- SOP 001. Laboratory Workflow and Documentation Procedures
- SOP 002. Balance Calibration
- SOP 003. Pipette Calibration
- SOP 004. Processing Urine Samples for Untargeted Metabolomics

Documentation is more than just keeping  
a lab notebook!





# The PDCA Cycle

<https://www.allaboutlean.com/pdca/>



# Case Studies QA Workshop

Annie Evans

# “Keep it real”

“Approval of a data set was delayed significantly because a single person used the wrong date format when they signed off on a final report... i.e. 5/24/18 rather than May 24, 2018. Everyone had to be tracked down again and resign the final report because there was an SOP which stated ‘use this date format’.”

# “Keep it real”

This highlights the need for ...

- Distinguishing between what is and is not important
- Keeping it simple!

When starting out with a Quality System, keep it within reason and grow as people become “used to” the process and mind set.

“Sure, I can pipet”

“We noticed a small and seemingly sporadic increase in variability of technical replicates in our QC review of data. We did a little digging and it turns out all the cases of increased variability were linked to a specific person preparing the samples. Turns out the individual was never properly trained on proper use of a pipet and so was pulling varying volumes when doing volume transfers”

“Sure, I can pipet”

This highlights the need for ...

- Training and continued proficiency testing
- Employee training logs

Also shows the power of good documentation as the lab was able to track the error back to a single person

## “Fear the wrath”

“Analyst saw QA department as predatory.... ‘they are trying to get me in trouble by reporting my errors’. This analyst never documented their deviations from SOP, because they were afraid of being ‘dinged’.”

# “Fear the wrath”

There is a helpful mindset with regard to QA and a not so helpful mindset.

- Not so helpful: THEY ARE OUT TO GET ME.. IT WILL GO ON MY RECORD...
- Not so helpful: YOU WILL BE COMPLIANT!!!!



# “Fear the wrath”

This highlights the need for ...

- A helpful mindset towards QA
- The need to collaborate with PIs on developing SOPs

QA may find things that will improve the process and help out with quality of data, but it is not intended to get in the way of getting science done

# “Consistent but wrong”

“A team of scientists was using a standardized excel macro to do a series of calculations to determine concentration. This was great because it built consistency into the workflow and made sure the team was using the same calculations in all their different projects....but no one had double checked the macro and the macro was doing some of the calculations incorrectly.”

# “Consistent but wrong”

This highlights the need for ...


- Properly developed SOPs that include system checks and acceptance criteria
- QA review from an outside perspective

A great mindset for a quality system, in general, is to think of it as a second set of eyes.

“You get what you pay for...”

“One of the first few studies we did was a large scale study of ~700 human samples trying to identify biomarker from cigarette smoking. Samples were sent to a metabolomics core facility, which hadn't run many large scale studies. We found several issues from the data...”

11369\_3



# Issues with the Retention Time Shift

Switching columns induced a retention time shift between samples

“You get what you pay for...”

This highlights the need for ...

- Properly developed SOPs that include system checks and acceptance criteria
- Proper training of instrument operators

At least good documentation allowed researcher to trace problem back to bad columns

The screenshot displays the ACQUITY UPLC Console interface. The left sidebar shows a tree view with 'Column Manager' selected and highlighted with a red box. The main area shows details for 'column ACQUITY UPLC@ BEH C18 1.7µm'. It includes usage counters for 'Injections Column 4' (467 injections), 'Total Samples on Column' (268 samples), and 'Total Sample Sets on Column' (26 sample sets). The environment section shows 'Maximum Pressure: 4/28/2011 10225 psi' and 'Maximum Temperature: 5/23/2011 60.3 °C'. The 'column history' table at the bottom lists various runs, with row 16 highlighted in red.


#	Date Started	Sample Set Name	User Name	System Name	Injections	Samples	Max psi	Max °C
11	4/27/2011 6:07 P	20110425	Administrato	INFORMATICS	87	45	8168	40.03
12	4/28/2011 12:32	20110425	Administrato	INFORMATICS	80	42	10226	40.04
13	4/29/2011 3:17 A	20110425	Administrato	INFORMATICS	8	5	8433	40.01
14	4/29/2011 4:45 A	20110425	Administrato	INFORMATICS	1	1	8437	40.01
15	4/29/2011 4:56 A	20110428	Administrato	INFORMATICS	3	3	8413	40.01
16	4/29/2011 5:30 A	20110428	Administrato	INFORMATICS	26	26	8487	40.01
17	5/3/2011 4:06 PM	20110502	Administrato	INFORMATICS	1	1	8744	40.01
18	5/3/2011 4:22 PM	20110502	Administrato	INFORMATICS	1	1	8694	40.01
19	5/4/2011 11:17 A	20110502	Administrato	INFORMATICS	1	1	8498	40.18
20	5/4/2011 11:28 A	20110502	Administrato	INFORMATICS	1	1	8378	40.08
21	5/4/2011 11:39 A	20110502	Administrato	INFORMATICS	21	11	8780	40.03
22	5/4/2011 4:14 PM	20110502	Administrato	INFORMATICS	1	1	8053	40.01
23	5/5/2011 11:42 A	20110502	Administrato	INFORMATICS	1	1	8662	40.3
24	5/5/2011 12:05 P	20110502	Administrato	INFORMATICS	1	1	8055	40.05
25	5/5/2011 12:17 P	20110502	Administrato	INFORMATICS	21	11	8563	40.03
26	5/23/2011 5:44 P	20110523	Administrato	INFORMATICS	1	1	6593	60.11

10980\_2

09-May-2011  
18:15:25

20110427PS\_ALL\_POS\_018


1: TOF MS ES+  
321.13 0.50Da  
67.6



10980\_1

20110427PS\_ALL\_POS\_661


1: TOF MS ES+  
321.1 0.50Da  
2.68e3



10980\_1\_2freethaw-cycle

20110427PS\_ALL\_POS\_637


1: TOF MS ES+  
321.13 0.50Da  
99.7



10980\_2

20110427PS\_ALL\_POS\_641

1: TOF MS ES+  
321.13 0.50Da  
85.1



10980\_3

# Issues in the Freeze-Thaw cycle

First, third and fourth traces are from same sample, freshly made and run. The second is after two freeze-thaw cycles. The ion 321.13 was not found in any of the fresh made samples, but appears after two freeze-thaw cycle.

“You get what you pay for...”

This highlights the need for ...

- Properly developed SOPs related to sample handling and storage
- Communication of quality procedures to clients upfront

Overall, the original data, which cost the client \$15k, was unusable, and the whole project had to be re-run (at a discounted price...)